

How to use your nebulizer with facemask



Neighborhood
Health Plan™



1. Wash your hands




2. Put the nebulizer machine on a flat, steady surface




3. Plug into outlet



4. Connect one end of the tubing to the nebulizer



5. Add medicine **and** saline
OR
Add premixed medication alone →



6. Attach facemask to the medication holder



7. Connect the other end of the tubing to the medication holder



8. Turn on the nebulizer machine



9. Place facemask over nose and mouth. Take **slow deep** breaths through your **MOUTH**



WRONG! Do NOT hold in front of face



10. When finished, wash parts in soapy water and rinse



11. Allow parts to air dry